
Winkelstabile Distale Humerusplatte

Angle-stable distal
humerus plate

Titan / Titanium

Variabel winkelstabil / Variable angle-stable

Königsee Implantate

www.koenigsee-implantate.de

Einleitung

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Fractures of the distal humerus have complex treatment requirements, especially when the articular surfaces are implicated. The angle-stable Königsee plates, which are positioned at 90° to one another, result in high primary stability of the osteosynthesis at the distal humerus. Consequently, functional after-treatment is generally possible, even in the case of osteopenic bone, and the risk of the injured joint stiffening is reduced.

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The implants are indicated for the surgical treatment of the following extra- and intraarticular fractures of the distal humerus according to the AO classification:

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- *Articular simple 13 C1.1-3*
 - *Articular simple, metaphyseal complex 13 C2.1-3*
 - *Articular multiple, metaphyseal simple – complex 13 C3.1-3*
- as well as pseudarthrosis and double fractures of the distal humerus*

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Fractures of type 13 A1.1-3 are only treated with plates in exceptional cases. 13 B3.1-3 fractures are stabilized with alternative implants.

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- *Bone and soft tissue infection in the region to be treated.*
- *Fractures in children and adolescents in whom epiphyseal fusion has not yet taken place.*

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General information about the treatment of distal humerus fractures

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The treatment of distal humerus fractures should be performed by experienced surgeons. The clinical examination includes the evaluation and documentation of the local skin and soft tissue situation, the recording of the sensorimotor neurological findings and of the arterial and venous vascular status. The physical examination is carried out with a focus on further potential injuries, particularly additional injuries of the injured upper extremity (caution > multiple fractures).

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Präoperativ-axial
Axial X-rays preoperative

Pre-operative imaging includes conventional X-rays of the elbow joint in 2 planes, spot-film radiography of the region if required, the demonstration of the adjacent bones and joints in 2 planes and a thin-layer CT of the elbow joint with reconstruction.

In the case of complex fractures within the context of multiple injuries, consideration should be given to a two-stage procedure with primary fracture immobilization by means of a cast or external fixation.

The internal osteosynthesis is performed once the soft tissue and the general state of health have been stabilized.

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The pre-shaped anatomical plate design increases the security of the application and reduces the time required for the operation.

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Two plates positioned at 90° to one another result in high primary stability of the osteosynthesis.

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The use of just one screw dimension reduces stock-keeping costs, increases the safety of the treatment and simplifies the reprocessing of the autoclavable container for implants and instruments.

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The ulnar plate screws can be positioned with variable angle stability. This simplifies the operation and reduces the risk of screws colliding.

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The articular block is fixed as far as the ulnar end of the trochlea of the humerus with a screw in the extension of the radial plate.

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
Screws in the radial and the ulnar plate work in opposite directions to address the main fragments of the epiphysis/metaphysis and ensure high primary stability.

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Longer plates are available on request for double fractures of the humerus.

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
The plates are anatomically pre-shaped and are available in various lengths. The radial plate is positioned dorsal-radial, and the ulnar plate ulnar-medial. Plate versions are available which are adapted for the relevant treatment side. The right plate is marked "R" on the top and the left plate is marked "L".

The plates must be combined only with designated angle-stable screws by Königsee Implantate.

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In the case of bilateral plate treatment, one of the two plates selected should be 2 holes shorter in order to reduce stress in the proximal bone and prevent the risk of adjacent fractures.

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The plates are fixed with angle-stable 3.5 mm screws with a 6 mm screw head. The screw holes of the radial plate are of an angle-stable design in order to prevent screws from colliding. In the ulnar plate, the screws can be inserted with variable angle stability – deviating from the standard direction of the screws by up to 15°. The screw hole position of both plates is designed in such a way in the humerus metaphysis / diaphysis that, when the plates are correctly positioned, the screws can be put in place without colliding. The thickness of the plate increases proximal to the humerus metaphysis. The undercut of the radial plate reduces the supporting surface at the bone and makes it easier to correct the plate shape. The implant design reflects the average anatomy. Bending instruments can be used to alter the shape of the plate. The bending process weakens the material and should be performed with caution and, wherever possible, once only. The direction of the screw in the region of the humerus diaphysis diverges alternately. This minimizes notch tension in the bone and the risk of fractures induced by the drilled hole and secondary implant migration.

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Positioning is possible in either a dorsal position or a stable lateral position over an upper arm bench with free mobility of the elbow joint.

The region to be treated must be fully demonstrable intraoperatively in 2 planes with the image converter. Perioperative antibiotics prophylaxis is obligatory. Partial deprivation of the upper arm blood supply is optional.

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In the case of extraarticular fractures, the approach with preservation of the extension apparatus should be considered. In the case of articular fractures, the V-shaped = Chevron osteotomy of the olecranon is recommended, together with dorsal capsule detachment and ulnar and radial soft tissue loosening of the triceps tendon in a proximal direction with extensive longitudinal, subtle exposure of the ulnar nerve. If the plates are to be fixed to the humerus diaphysis in a very proximal position, the radial nerve is subtly exposed.

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Firstly, in the case of articular fractures, the trochlea of the capitulum and the trochlea of the humerus are anatomically reduced. K-wires or transversely positioned screws in the exact center of the trochlea of the humerus – in order to prevent subsequent collision with plate screws – fix the articular block. Axial compression must be avoided at all costs. If there are osseous defects, consideration should be given to building up the bone with autologous bone or bone substitute materials. Small osseous and osteochondral fragments are fixed temporarily with fine K-wires. In the next step, the articular block is fixed to the metaphysis / diaphysis with K-wires with the correct anatomical axial alignment and rotation.

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With regard to the surgical approach, we would point to the relevant literature.

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
We would point to the surgical steps outlined in surgical teaching. Please follow the generally applicable rules for osteosynthesis.

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
Depending on the fracture situation, the definitive unilateral or bilateral plate osteosynthesis is now performed. Select a plate of a suitable length and of the correct side. A minimum of 6 in case of long plates 8 cortical screw contacts should be inserted into the humerus metaphysis/ diaphysis on the other side of the proximal fracture line. We recommend bicortical screw insertion in the humerus diaphysis, corresponding to 2 cortical screw contacts. In the case of bilateral plate fixation, start by fitting the radial plate.

The plate osteosynthesis itself is described in the following chapters "Radial plate" and "Ulnar plate". Finally, smaller osseous and osteochondral fragments are fixed definitively at the articular block using biodegradable pins, suture anchors or small flat-head screws. The implants should lie at or below the cartilage level. Following the temporary refixing of the osteotomized olecranon, the play of the joint, the joint axes and the trouble-free motion sequence of the joint are checked clinically and radiologically.

A radiological check of the findings with an image converter in at least 2 planes is obligatory.

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
The plates are anatomically pre-shaped and are fixed using angle-stable 3.5 mm screws with a 6 mm screw head. Please note that there is a right version marked "R" on the top of the plate and a left version marked "L", available in various lengths. The proximally positioned section of the plate has undercuts which make the reshaping of the plate easier. The screw direction in this area of the plate diverges.


The screws in the distal section of the radial plate address the capitulum of the humerus with three screws. A screw positioned in the radial plate extension addresses the distal, ulnar end of the trochlea of the humerus.

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Articular fracture elements can be fixed by means of transversely inserted screws or wires in the area marked with a circle without any implant collision.

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

Select a suitable plate and check that the side – right or left – is correct and of a sufficient length for the fracture morphology. Please note that, for bilateral treatment, one of the two plates selected must be 2 holes shorter. Screw one angle-stable drill guide into the radial extension and another one into one of the two distal plate holes. Now push the plate onto the humerus from a distal direction. The distal end of the plate must lie securely proximal to the bone-cartilage junction, in order to prevent any striking of the head of the radius. The distance between the distal end of the humerus and the end of the radial plate is between 5 and 10 mm. A recess in the underside of the distal end of the plate hooks into the end of the bone when the plate is correctly positioned. Make sure that the alignment to

the humerus shaft is correct. Now push the plate in an ulnar direction until the radial plate extension fits closely against the radial epicondyle. The implant must not encroach into the olecranon fossa.

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The plate shape reflects the average anatomy. If there is insufficient contact between the plate and the bone, you can use conventional instruments to carefully reshape the plate once; the bending area is shown with α . As a general rule, the bending movement is performed towards the underside of the plate. According to the anatomical preliminary examinations, this may be necessary for humeri with small dimensions. Please note that bending processes weaken the plate material with the associated risk that the implant may break.

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Hold a K-wire parallel to the drill guide inserted at the distal end of the plate and check whether the capitulum of the humerus can be successfully addressed with screws. Now place the K-wire parallel to the drill guide inserted in the plate extension and check whether the ulnar part of the trochlea of the humerus can be addressed with a screw.

Prior to the drilling process, make sure that the plates are aligned correctly to the humerus shaft. Fix the plate manually with OP assistance, reduction forceps or the Königsee REP-RET system.

First insert a screw into one of the two distal plate holes. Predrill through the drill guide using a 2.5 mm drill. In the case of osteoporotic bone, the cortical bone close to the plate can be prepared with a 2.5 mm drill and the rest of the screw channel can be prepared with a long 2.0 mm drill. The insertion of the screw compresses the bone around the screw and increases the primary stability between the screw and the bone. Determine the length of the screw, cut the thread and screw in a 3.5 mm cortical screw of this length. There is no need to cut the thread if a self-tapping screw is used. Interlock the screw head thread with the plate hole thread without applying excessive force in order to prevent cold welding between the screw head and the plate hole. In the same way, insert screws in additional screw holes in the distal part of the plate and check the position of these using anterior-posterior, lateral and, if necessary, oblique image converter control. In the proximal part of the plate, screws are inserted bicortically.



The screw direction diverges alternately. Proximal to the fracture, 3 screws which grip bicortically should be inserted for short plates and 4 for long plates.

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The subsequent procedure is the same as for the radial plate. The screw direction in the proximal part of the plate which fits closely against the humerus diaphysis diverges in order to reduce notch stress in the bone and increase the pull-out strength of the plate. Select a plate of sufficient length for the fracture morphology and please note that one of the two plates should be at least 2 holes shorter. This decreases the load on the bone and risk of secondary fractures. Check the fit against the bone and, if necessary, use conventional bending instruments to bend the plate. Insert screws in the plate holes in the distal part of the plate using angle-stable drill guides. Hold a K-wire in the flight of the drill guide and check whether the capitulum of the humerus can be addressed collision-free, as required, with the screws already inserted. If necessary, use the variable angle stability of the plate hole. The subsequent procedure is the same as for the radial plate.

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Once the osteosynthesis has been completed, the free play of the joint is checked clinically and radiologically. Impingement between the osteosynthesis material and the osseous structures should be prevented. Dissected nerves are positioned loosely in situ and the position in relation to the inserted material is documented in the OP report. Refix the osteotomized olecranon. The result is checked radiologically with the image converter and documented. Following the insertion of drains, soft tissue and skin are adapted over the implants without tension. Sensitivity, motoricity and the arterial and venous vascular status are checked postoperatively and documented.

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The osteosynthesis stability achieved and patient compliance limit the after-treatment and consideration should be given to temporary immobilization with a splint or orthosis. Physiotherapeutic after-treatment is introduced at an early stage. The removal of material is optional and should only take place once the fracture has healed securely and, at the earliest, 18 to 24 months following the osteosynthesis.

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
19.550.00

1

Perforated autoclavable container with inset for instruments and implants

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6.031.16

5

Kirschner wire with trocar point and round end, diameter 1.6 mm, length 150 mm

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2.904.06

1

Twist drill for quick coupler diameter 2.5 mm, length 125 mm

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2.953.60

1

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Screw gauge with clasp, for screws diameter 3.5 - 4.0 mm, MR 60 mm

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2.940.25SA

1

Hexagonal screwdriver, conical with handle,

for screws diameter 2.7 - 4.0 mm, wrench width 2.5

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*No part of this publication may be reproduced,
stored in a retrieval system or transmitted in any
form or by any means electronic,
mechanical, photocopying, recording or other-
wise without the prior written permission of the
publisher.*

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© Copyright Königsee Implantate GmbH
Königsee Implantate GmbH