Winkelstabile Distale Radiusplatten

Angle-stable distal radius plates

Titan / Titanium

Minifragment / Mini fragment

Königsee Implantate

www.koenigsee-implantate.de

Einleitung

Introduction

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90% of distal radius fractures are extension fractures. The difficulties of the dorsal approach and the added operative trauma of extensor tendon irritation and chiselling away Lister's tubercle have led in recent years to a change in the management of these fractures with a trend towards palmar angle-stable plate internal fixation.

However, from the biomechanical aspect, a plate applied dorsally undoubtedly offers considerable advantages.

The angle-stable radial plate for dorsal treatment was developed despite this conflict.

Because of the Y shape, Lister's tubercle with the extensor pollicis longus tendon running around it is left intact.

The 5 distal 2-mm screws are secured in the plate at a fixed angle and support the joint surface through their ordered anatomical arrangement. Even small fragments such as the typical dorso-ulnar border fragment can be fixed securely in this way.

The longitudinal hole lies in the proximal shaft area of the plate. After fine repositioning of the plate, the correct axis of the plate in the shaft region is guaranteed. The mechanical strength of the two distal shaft screws on which the greatest forces act is much higher as these are screwed in only once.

The design of the dorsal angle-stable radial plate ensures easy and secure treatment of especially dorsally tilted distal radius fractures, and the additional operative trauma is much lower compared to the previous procedure.

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The 2 mm palmar radius plate is an implant precisely designed for the anatomical conditions of the palmar distal end of the radius. The six distal 2 mm screws are fastened to the plate with angular stability and support the joint area in an anatomically aligned gallery arrangement. Small fragments supporting the joint area can also be securely fixed in this way.

The slot is in the proximal shaft region of the plate. It ensures that the plate is positioned correctly in the shaft region after the fine repositioning of the plate. The mechanical strength of the two distal shaft screws, which are subject to the greatest forces, is significantly higher, because they are only rotated once.

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

*Extension fractures of distal radius
comminuted fractures with joint involvement.*

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

- *Smith fractures, reversed Barton fractures (23 – B3)*
- *Distal radius extension fractures with and without involvement of the joint (23 – A3 / C1 / C2 / C3, with preference for palmar access)*

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*Supine position, hand table, bloodless field,
extension using traction extension sleeves on
the first and second digits (3-4 kg).*

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


*Exposure of the extensor retinaculum.
Incision of the 4th radial compartment of the tendons of the common extensor digitorum.*

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Dissect close to bone, elevating the extensor digitorum communis tendon in the ulnar direction.

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
Incision of the 2nd ulnar compartment of the tendons of the extensor carpi radialis brevis et longus.

Dissect close to bone, retracting the tendons in the radial direction; the 3rd tendon compartment with Lister's tubercle and the tendon of extensor pollicis longus are left intact.

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

Fine reduction of the joint surface and dorsal fragments with periosteal elevator, impactor and finger pressure.


Grafting with cancellous bone or bone substitute is rarely necessary with the angle-stable plate.

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*Underrunning the tendon of extensor pollicis longus, applying the plate and filling the middle of the longitudinal hole with diameter 3.5 mm cortical screw.
Fine repositioning of the plate, filling of the central shaft.*

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
Screw on the 2 mm drill guide, drill, remove the drill guide, measure the length and screw in the distal 2 mm screws maintaining the angle.

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
*Fill the distal shaft hole.
Carry out a final check under image intensifier
control and record the internal fixation
with printer.*

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*Suture of retinaculum, tourniquet release,
vacuum drain, wound closure in layers.*

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Postoperative brief immobilisation in a plaster splint, immediate exercises, ergotherapy etc. in the splint.

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
*Supine position, hand table, bloodless field,
extension using traction extension sleeves on
the first and second digits (3-4 kg).
Sterile cover. 7 cm long skin incision radial
from the flexor carpi radialis tendon.*

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

*Preparation between flexor carpi
radialis tendon and radial artery.
L-shaped removal of the pronator quadratus.
Exposure of fracture.*

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
Fine repositioning of the fracture. Fit plate and insert diameter 3.5 mm cortical screw in centre of slot. Fine correction of plate position, other shaft holes covered.

When restoring extension fractures the distal plate holes can also be first fixed with 2 mm screws for angular stability with a plate raised from the shaft. Then the repositioning is carried out by raising the plate at the radius shaft.

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Screwing in the 2 mm drill sleeve, pilot drilling, removal of the drill sleeve, determining the length, screwing in the 2 mm screws for angular stability. They must not protrude above the dorsal cortex.

As many screw holes as possible should be used depending on the fracture with its fragments. Final check under BW and documentation of the osteosynthesis with printer.

Open the blood block, Redon drainage, close wound.

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Cortical screws, diameter 2.0 mm, with conical head thread, fully threaded, self-tapping

Bestell-Nr. Titan <i>Code N° titanium</i>	Länge <i>Length</i>	Anzahl im Set <i>Quantity in set</i>
3.114.16	16 mm	5
3.114.18	18 mm	5
3.114.20	20 mm	5
3.114.22	22 mm	5
3.114.24	24 mm	5
3.114.26	26 mm	5
3.114.28	28 mm	5
3.114.30	30 mm	5

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Cortical screws, diameter 3.5 mm, fully threaded, self-tapping

3.132.14	14 mm	5
3.132.16	16 mm	5
3.132.18	18 mm	5
3.132.20	20 mm	5
3.132.22	22 mm	5
3.132.24	24 mm	5

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Angle-stable distal palmar radius plate for screws diameter 2.0/3.5 mm

5.898.63	rechts / right	50 mm	6/3	1
5.898.65	rechts / right	67 mm	6/5	-
5.899.63	links / left	50 mm	6/3	1
5.899.65	links / left	67 mm	6/5	-

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Angle-stable distal dorsal radius plate for screws diameter 2.0/3.5 mm

5.892.53	rechts / right	63 mm	5/3	1
5.893.53	links / left	63 mm	5/3	1

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Angle-stable distal dorsal radius plate, narrow for screws diameter 2.0/3.5 mm

5.892.43	rechts / right	63 mm	4/3	-
5.893.43	links / left	63 mm	4/3	-

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Kirschner wire with trocar point and round end, diameter 1.4 mm, length 150 mm

6.031.16

5

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D 1,5 mm, Länge 105 mm

Twist drill for quick coupler diameter 1.5 mm, length 105 mm

2.904.20

1

D 2,5 mm, Länge 105 mm

Twist drill for quick coupler diameter 2.5 mm, length 105 mm

2.904.05

1

Messinstrument

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Screw gauge for screws diameter 2.0 mm, MR 44 mm

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Hexagonal screwdriver, conical with handle,

for screws diameter 2.7 - 4.0 mm, wrench width 2.5

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*No part of this publication may be reproduced,
stored in a retrieval system or transmitted in any
form or by any means electronic,
mechanical, photocopying, recording or other-
wise without the prior written permission of the
publisher.*

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© Copyright Königsee Implantate GmbH
Königsee Implantate GmbH