
Variabel winkelstabile Olecranonplatte

Variable angle-stable
Olecranon plate

Titan / Titanium

Variabel winkelstabil / Variable angle-stable

Königsee Implantate

www.koenigsee-implantate.de

Introduction

Complex fractures of the proximal ulna frequently pose a considerable challenge for the trauma surgeon.

Following the evaluation of any secondary injuries, the aim should be the anatomical reconstruction of the joint structure, together with the secure anchorage of the extensor system. However, even though the complex proximal ulnar fracture is a relatively rare fracture, the proportion of osteoporotic cases is rising. When there is more osteoporotic bone, osteosynthesis with conventional implants is associated with a higher failure rate.

The aim of the new, variable, angle-stable olecranon plate is to increase the primary stability that can be achieved and to make surgery easier thanks to the anatomically pre-shaped design of the plate.

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As a general rule, the variable, angle-stable olecranon plate can be used to treat both simple and complex extra- or intraarticular fractures of the proximal ulna (Morrey Type I A to Type III A;B). Further indications are pseudarthroses or osteotomies of the proximal ulna.

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Type I: non-dislocated fractures (dislocation < 2mm):
I A: single fragment
I B: several fragments
Type II: dislocated but stable fracture (dislocation > 2mm, but intact collateral ligaments):
II A: single fragment
II B: several fragments
Type III: dislocated and non-stable fracture
III A: single fragment
III B: several fragments with possible fracture of the radial head

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The contraindications are essentially the same as those for conventional plate osteosynthesis, such as acute infections, childhood fractures, non-stable soft tissue conditions or insufficient soft tissue coverage in very thin patients.

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As well as taking the current case history and the findings of a clinical examination, X-rays in two planes are generally sufficient for the pre-operative diagnosis. In the case of complex fractures, a CT with accompanying reconstruction may be helpful.

Any secondary injuries, such as radial head fractures or Monteggia fractures must be identified and treated additionally as required.

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The proximal ulnar region has just a small amount of soft tissue coverage. The olecranon has a tendency towards dislocation due to the pull exerted by the triceps tendon. The trochlear articular surfaces are often implicated in the fracture. An implant in the proximal ulnar region therefore needs to be able to meet these special requirements.

The new, variable, angle-stable plate corresponds to a small-fragment DC plate in the ulnar shaft region and has both variable angle-stable holes and DC holes. Since the middle section of the ulna often has a slightly arched shape, the plate can be adapted to the ulna by means of two reconstruction profiles.

The plate becomes broader as it approaches the trochlear region, is bent slightly and is thinner here, in order to prevent it from irritating the soft tissues in this area.

A total of four variable angle-stable holes have been provided in a slightly offset arrangement.

A further two 1.5 mm holes have been positioned in this section, either for the temporary fixation of the plate with Kirschner wires or for the mounting of suture material for the trans osseous fixation of smaller fragments.

The plate section for the olecranon is in the center. The proximal plate section is bent in a fixed position and contains two variable angle-stable holes. This special design means that, if desired, a long screw can be inserted that extends into the shaft.

Two design versions of the bent section are available: one version has two claws at the head for additional stability, while the other version has two eyelets at the head, meaning that two interlacing sutures through the triceps tendon using non-absorbable fiber-reinforced suture material can be securely anchored to the plate. This enables the tensile forces that occur to be safely neutralized.

The plate is available in various lengths. Since the plate can be adapted, there is no need for a separate left and right version, which means that storage costs are reduced.

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
Pre-operative: 91-year-old female patient, fall onto the right arm

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Treated surgically with an angle-stable olecranon plate, design version with holes for reinforcement of the triceps tendon with suture material.

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85-year-old female patient, fall onto the right arm, additional multi-fragment fracture of the humeral head.

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

Treated surgically with an angle-stable olecranon plate, design version with claws.

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
The patient is placed in a supine position with the arm positioned on the upper body or, alternatively, in a prone position with the forearm hanging down.

The fracture is exposed by means of a dorsal approach along the radial ulnar edge with radial circumscision of the tip of the olecranon. The incision can be extended in a distal direction – as required by the extent of the fracture. The anconeus muscle is firstly stripped away on the radial ulnar side, together with the radial joint parts, if required.

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Following the appropriate reduction of the individual fragments, these are then fixed with Kirschner wires as temporarily as possible. Smaller fragments with attached ligament parts can also be reinforced with appropriate suture material and indirect reduction can therefore be achieved.

Impressions of the articular surface are taken and cancellous bone lining can be performed if required. Wherever possible, the Kirschner wires should be positioned in such a way that they do not interfere with the plate installation. Since the ulna is generally slightly curved, the plate can be bent slightly in the region of the reconstruction profile if necessary.

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
If the surgeon wants to position the part of the plate that rises towards the tip of the olecranon directly against the bone, the triceps tendon must firstly be cut open lengthwise and cut away from the tip of the olecranon in this dorsal part of the olecranon.

A small space between the plate and the bone in the proximal ulnar region is acceptable.

With the two available versions of plate design - with and without claws - the tensile forces of the triceps tendon are transmitted to the plate in different ways: either via the hook-shaped ends of the plate or by means of anchor sutures in the triceps tendon which are pulled through the suture fixation holes in the alternative plate design.

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
The triceps tendon is now interlaced in this region with PDS II sutures or PDS tapes which are then pulled through the two suture holes in the ascendant part of the plate. The plate is then put in position and is first fixed to the distal ulna by the moderate tightening of the screw in the longitudinal hole. Following appropriate fluoroscopy, the position of the plate is corrected if required. The screw in the longitudinal hole of the plate is then tightened securely.

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A screw is then inserted from a dorsal position into the fracture fragment supporting the triceps tendon; the screw is inserted at a slight upwards angle and grips the opposite cortical layer in the region of the base of the coronoid process. A minimum of two - and ideally three - cortical screws

should be placed distally to the main fracture. Either standard or angle-stable screws can be used here, depending on the osteosynthesis desired. The screws should be inserted into the angle-stable screw holes in the actual joint region in as divergent a manner as possible in order to prevent screw irritation.

The sutures looped through the triceps tendon insertion point act like a protective tensioning band for the triceps tendon and these are now knotted at the plate.

The final step of the osteosynthesis procedure is intraoperative fluoroscopy in order to ensure that there is no intraarticular osteosynthesis material left. A functional check for existing exercise stability should also be performed. The joint capsule is closed with absorbable suture material and this is followed by the possible insertion of a Redon drain and the appropriate closure of the wound.

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During surgery, a protective upper arm splint is applied. After-care depends on the osteosynthesis stability achieved, the bone quality and patient compliance.

Immobilization until the wound has healed is the usual procedure, with guided exercises while wearing the splint. However, a more prolonged period of immobilization may also be required depending on the level of stability achieved.

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
2

Kirschnerdraht Trokar, Ende rund, D 1,6 mm,
 Länge 150 mm

Kirschner wire with trocar point and round end, diameter 1.6 mm, length 150 mm

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D 2,5 mm, Länge 125 mm

Twist drill for quick coupler diameter 2.5 mm, length 125 mm

2.904.06

1

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Screw gauge with clasp, for screws diameter 3.5 - 4.0 mm, MR 60 mm

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Hexagonal screwdriver, conical with handle,

for screws diameter 2.7 mm - 4.0 mm, wrench width 2.5

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*No part of this publication may be reproduced,
stored in a retrieval system or transmitted in any
form or by any means electronic,
mechanical, photocopying, recording or other-
wise without the prior written permission of the
publisher.*

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© Copyright Königsee Implantate GmbH