
Variabel winkelstabile **Patella-Platte**

Variable angle-stable
Patella-plate

Titan / Titanium

Variabel winkelstabil / Variable angle-stable

Königsee Implantate

www.koenigsee-implantate.de

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
In a number of biomechanical studies, the angle-stable, bilateral patellar plate shows a significantly higher level of stability and rigidity in transverse patellar fractures in comparison with traditional (modified) tension-band osteosynthesis with K-wires and the cannulated screw osteosynthesis with anterior tension band. In the case of tension-band osteosynthesis, the use of angle-stable, bilateral patellar plate osteosynthesis can - due to its high level of rigidity - eliminate the risk of the commonly observed secondary dislocation of the fragments as a result of soft tissue atrophy and the degradation or cutting of the stainless steel wire during post-operative treatment, and the primary reduction result can be better retained. Broken cerclage wires and K-wire migration, which make it necessary to remove material prematurely in 10-65 % of cases, can also be avoided through the use of bilateral, angle-stable plate osteosynthesis.

In the case of patella fractures, angle-stable bilateral plate osteosynthesis requires the surgeon to rethink his approach due to the anatomy of the patella and the biomechanical functional principle, which differs from that of tension band osteosynthesis. Bilateral angle-stable plate osteosynthesis is based on the biomechanical principle of internal fixation. This means that, for biomechanical reasons, a radiologically detectable distance between the plate and the patella is desirable rather than something to be avoided, as this ensures that the soft tissue surrounding the patella (quadriceps tendon, patellar tendon, retinacula) is protected.

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- Multi-fragment patellar fractures
(if there is a sufficiently large proximal pole fracture for an angle-stable screw to be safely introduced in each case)
- distal patellar pole fractures
(if there is a sufficiently large distal pole fragment, which is safely covered by the distal portion of the plate)
- Comminuted patellar fractures
without the possibility of secure anchorage for screws
- Proximal patellar pole fractures
- Osseous detachment of the patellar tendon

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

The operation is carried out in the dorsal position and the contralateral leg is slightly lowered to enable an overlap-free intraoperative fluoroscopy in the lateral ray path. Positioning the affected knee joint in the extended position facilitates digital monitoring of the posterior patella during the reduction and during the drilling process, on account of the lower contact pressure on the patella. In the case of obese patients, slight flexion of the knee joint by supporting the back of the knee with a small rolled-up towel can help to prevent obstruction caused by skin and subcutaneous tissue. A thigh tourniquet may be applied.

Zugang Access

Zugang zur Patella
Access to the patella

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
In principle, it is possible to access the patella both longitudinally and transversely. Although transverse access gives a slightly better cosmetic result, longitudinal access to the patella offers several advantages: Less obstruction is caused by skin and subcutaneous tissue – especially in the case of adipose patients, access can be extended easily if necessary and the access does not compromise subsequent implantation of a total knee endoprosthesis, particularly in critical soft tissue conditions.

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

Firstly, the longitudinal incision of the skin and subcutis is made directly over the patella, starting 2 cm above the upper patellar pole to a length of approx. 10 cm in a distal direction as far as the level of the distal patellar pole. This is followed by deeper dissection as far as the prepatellar bursa, which, if applicable, is resected in the case of a rupture. Exposure of the medial and lateral supporting structures of the patella (retinacula of the patella).

To ensure that the plate is positioned securely, it is recommended that the medial and lateral edge of the patella is exposed with the retinaculum to a width of two to three centimetres.

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If the retinaculum is not ruptured, the intraarticular fracture hematoma can be drained by means of a longitudinal, lateral arthrotomy, a finger's breadth lateral to the edge of the patella, and the retropatellar surface can be palpated to check the reduction. The reduction of the patellar fracture is monitored digitally and the reduction result is secured with the aid of reduction forceps. In the case of a transverse patellar fracture, it is advisable to position the reduction forceps at 90 degrees to the fracture line, centrally above the patella, from the upper patellar pole to the tip of the patella. The application of the bilateral, variable, angle-stable patellar plate is thus not hindered by the positioning of the reduction forceps.

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

Manual intraoperative bending of the arm

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When bending the plate, remember that repeatedly bending the delicate titanium plate forwards and backwards leads to material damage, so the plate must be adjusted by bending it gradually and only in one direction. One should avoid bending the plate back.

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
Long-nose forceps are used to bend the small arms close to the base, so that they enclose the distal patellar pole in a dorsal direction, until a slight hook configuration is achieved.

Both plate shafts can be adjusted by bending between the plate holes, to match the patella's radius, using two screw drill sleeves or two forceps.

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After positioning two longitudinal stab incisions medial and lateral to the patellar tendon at the lower patellar pole, hook the arms medially and laterally next to the patellar tendon while enclosing the distal patellar pole, and rest the patellar plate on the patella.

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Stab incisions to position the patellar plate

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Here, the two plate shafts should come to rest at an angle of 70-90° to the frontal plane at the lateral and medial edge of the patella and both caudal arms should securely enclose the distal patellar pole while leaving out the patellar tendon.

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Gradual remodeling of the patellar plate may be necessary in order to adapt it to suit the anatomy. Under no circumstances should both plate shafts be bent 90° to the frontal plane, as the screws will then come to rest ventrally and will not grip sufficient osseous tissue. Although both arms come to rest on the interior of the patella so that they can grip the distal patellar pole, there is no reason to fear an intraarticular position or impairment of the hyaline articular cartilage, because in physiological terms, the distal patellar pole does not have any hyaline articular cartilage covering and does not articulate with the trochlea.

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
After the patellar plate has been fitted, either both plate shafts can be braced against one another with pointed reduction forceps via the eyelets lying at the cranial end of the plate, or in each case one plate shaft can be fixed to the patella.

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
Low pressure should be applied when closing the reduction forceps, so that dislocation of the plate over the convexly curved patella does not occur.

Depending on the position of the main fragment, first of all two screws of corresponding lengths are inserted into either the most proximal (in the case of a proximal main fragment) or the most distal (in the case of a distal main fragment) plate holes.

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
It is advisable here to create the first hole so that it is angle-stable via the threaded drill guide, provided that the position of the plate guarantees that the screw to be inserted crosses the fracture line and preferably passes through the entire patella.

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Pre-drilling the screw holes with the aid of an angle-stable or variable angle-stable drill guide.

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Otherwise, the hole should be made using the variable, angle-stable drill guide to ensure that the screw is in the correct position. It is advisable to create the hole using intraarticular digital monitoring of the lateral arthrotomy, in order to avoid drilling into the joint. After measuring the length, a screw of appropriate length is inserted. The screw head should be locked into the plate thread slowly and with a small amount of torque, as otherwise a loss of reduction may occur.

Subsequently, a screw of corresponding length is inserted into the most proximal or most distal plate hole lying opposite using the variable, angle-stable drill guide, and the screw should also cross the fracture line and pass through the entire patella. In order to avoid colliding with the previously inserted screw, the hole should be drilled in a different direction to that of the previously inserted screw using variable angular stability.

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After the holes have been drilled and the length has been measured, screws of corresponding length are then inserted into the remaining plate holes. Only those plate holes which are located directly at the level of the fracture do not have screws inserted into them. A short screw can optionally be inserted into the ventral plate hole above the distal patella, as long as the distal patellar fragment is large enough. When filling this plate hole, perforation of the subchondral lamella should be avoided, in order to prevent an intraarticular position of the screw. This screw usually does not have to be inserted. After a screw has been positioned respectively proximal and distal to the fracture in each plate shaft, the reduction forceps can be removed. All screws in the plate should be additionally tightened and anchored alternately and gradually, to avoid secondary dislocation. The final digital intraarticular palpation should include intraarticular, step-free reduction without dehiscence and avoid intra-articular screw malposition.

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
During fluoroscopy monitoring in both planes, sufficient screw length should be ensured and the screw should not be positioned intraarticularly, wherein the digital monitoring of the posterior patellar surface is often more significant than fluoroscopy monitoring on account of the shape of the patella. A dehiscence of several mm between the plate and the edge of the patella is normal and desirable.

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In order to avoid colliding with the previously inserted screws when drilling, the hole should not be too flat and the position of the plate at the edge of the patella should not be too lateral with plate holes lying directly opposite to one another. Due to the increasing thickness of the patella in a medial direction, the holes should be drilled - under digital monitoring - through the patella from the lateral side at an angle of approximately 10-15 degrees to the frontal plane, while the medial holes should be drilled at an angle of 0 degrees.

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
To ensure the stability of the plate osteosynthesis, it is sufficient for the two screws to cross the fracture line in a divergent or convergent manner and, if possible, to cross the patella over its entire diameter. All other screws can be inserted just a short distance and must not cross the fracture line.

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Any torn retinacula should be sutured immediately after repeatedly rinsing the wound, as this significantly increases the stability of the osteosynthesis. Following placement of intraarticular drainage with suction in the upper recessus, the closure of the lateral arthrotomy is carried out. Subcutaneous drainage with suction should be avoided as this could lead to pressure necrosis of the contused thin skin over the patella. After the layers of skin tissue have been closed, a sterile elasticated dressing is applied to the leg.

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Simple transverse patellar fractures with two large main fragments can receive early functional treatment with a partial loading of 20 kg over six weeks after the operation. Full loading may be permitted successively, depending on the results of the X-ray. Distal pole patellar fractures, multi-fragment patellar fractures and longitudinal patellar fractures should be treated on an individual basis and according to the fracture pattern. In these cases, the level of mobility can be limited (e.g. two weeks at 30° flexion, two weeks at 60° flexion, two weeks at 90° flexion), if necessary, for six weeks after the operation in an appropriate orthosis.

- S. Thelen, J. Schneppendahl, E. Jopen, C. Eichler, J. Koebke, E. Schonau, M. Hakimi, J. Windolf, M. Wild. Biomechanical cadaver testing of a fixed-angle plate in comparison to tension wiring and screw fixation in transverse patella fractures. *Injury* (2012), <http://dx.doi.org/10.1016/j.injury.2012.04.020>
- S. Thelen, J. Schneppendahl, R. Baumgartner, C. Eichler, M. Betsch, M. Hakimi, J. Windolf, M. Wild. Cyclic long-term loading of a bilateral fixed-angle plate in comparison with tension band wiring with K-wires or cannulated screws in transverse patella fractures. *Knee Surgery, Sports Traumatology, Arthroscopy* (2012)
- M. Wild, S. Thelen, P. Jungbluth, M. Betsch, D. Mietsch, J. Windolf, M. Hakimi. Fixed-Angle Plates in Patella Fractures – A Pilot Cadaver Study. *Eur Med Res* 2011, 16(1) : 41 - 46
- M. Wild, C. Eichler, S. Thelen, P. Jungbluh, J. Windolf, M. Hakimi. Fixed-angle plate osteosynthesis of the patella – An alternative to tension wiring? *Clin Biomech* (2010), 25(4) : 341 - 47.

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
Cortical screws, diameter 3,5 mm, with conical head thread, fully threaded, self-tapping

Bestell-Nr. Titan Code N° titanium	Länge Length	Anzahl im Set Quantity in set
3.136.12	12 mm	5
3.136.14	14 mm	5
3.136.16	16 mm	5
3.136.18	18 mm	5
3.136.20	20 mm	5
3.136.22	22 mm	5
3.136.24	24 mm	5
3.136.26	26 mm	5
3.136.28	28 mm	5
3.136.30	30 mm	5
3.136.32	32 mm	5
3.136.36	36 mm	5
3.136.40	40 mm	5

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Perforated autoclavable container with inset for instruments and implants

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Drill guide for small fragment, for variable angle-stable screwing, length 45 mm

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Drill guide for small fragment, for angle-stable screwing, length 45 mm

2.977.13

2

Kirschnerdraht Trokar, Ende rund,
D 1,0 mm, Länge 150 mm

Kirschner wire with trocar point and round end, diameter 1.0 mm, length 150 mm

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D 2,5 mm, Länge 125 mm

Twist drill for quick coupler diameter 2.5 mm, length 125 mm

2.904.06

1

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Screw gauge with clasp, for screws diameter 3.5 - 4.0 mm, MR 60 mm

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*Hexagonal screwdriver, conical with handle,
for screws diameter 2.7 mm - 4.0 mm, wrench width 2.5*

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*No part of this publication may be reproduced,
stored in a retrieval system or transmitted
in any form or by any means electronic,
mechanical, photocopying, recording or
otherwise without the prior written
permission of the publisher.*

SUPPLIED BY

ATHRODAX
HEALTHCARE INTERNATIONAL LTD

HAWTHORN BUSINESS PARK, DRYBROOK
GLOUCESTERSHIRE, GL17 9HP
Telephone: (01594) 544440 Fax: (01594) 545800
Email: sales@athrodax.co.uk www.athrodax.co.uk

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© Copyright Königsee Implantate GmbH
Königsee Implantate GmbH