Winkelstabile Distale Tibiaplatte

Kleinfragment

Angle-stable distal tibia plate
Small fragment

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The management of distal fractures of the tibia, particularly those involving the joint surfaces, is demanding. The implant presented here combines the advantages of a fixator with the principle of an angle-stable plate that allows internal fixation. The distal tibial plate is available in different lengths and is shaped. It is fixed to the tibia medially.

The implant is designed in small-fragment dimensions. The distal tab provides a fixation of the medial malleolus with a conventional small-fragment screw. At the distal end of the plate, there are angle-stable screw holes that allow a variation of up to 20° in the angle of the screws. The given radiation of the screws ensures that the screws can achieve angle-stable fixation of all the main fragments of the epiphysis and metaphysis without touching. There are two more screw holes on the anterior plate flange which are used for fragment fixation in the antero-posterior direction. The plate cross-section is reinforced proximal to the metaphysis. The undercut diminishes the contact surface on the bone and allows the plate shape to be corrected. The implant reproduces the average anatomy and the shape can be altered with the usual bending instruments. On the proximal limb of the plate, the surgeon can choose to occupy the screw holes with conventional or angle-stable small-fragment screws because of the special design (combi-hole). When inserted as angle-stable screws, the screws diverge alternately. This diminishes notching stresses in the bone and thus the risk of secondary implant migration.

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The implant was developed for the operative management of fractures of the intra-articular distal tibia and also extra-articular fractures of the distal tibia.

If a fracture running in the frontal plane is present, (additional) anterior plate stabilisation by a different implant is necessary.

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General remarks on the management of fractures of the distal tibia

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The management of distal tibial fractures, especially intra-articular fractures, is demanding and should be reserved for experienced surgeons. With complex fractures, a two-stage procedure is indicated with the primary treatment of the fibula, minimally invasive reduction of the joint block and fixator assembly. Internal fixation is performed only when the soft tissue has been stabilised. Preoperative imaging includes conventional radiographs of the ankle joint in 3 planes, of the lower leg and knee in 2 planes and a thin-layer CT if the intra-articular distal tibial shaft is involved. The perfusion of the soft tissues, particularly if there is a concomitant fibular fracture with the need for a second incision, and existing soft tissue damage by the trauma pose high demands in the choice of approach.

Immediate medial incisions should be avoided and the relevant literature should be referred to. Implant coverage must be ensured. An alternative is the coverage with a free flap within 48 hours after internal fixation. The patient is in supine position and is given perioperative antibiotic prophylaxis; use of a thigh tourniquet is optional.

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The shaped anatomical plate design increases the safety of use and reduces the operation time.

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Use in articular and extra-articular tibial fractures.

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Use of the small-fragment size alone reduces stock maintenance costs and facilitates tray preparation.

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The angle-stable course of the screws in the metaphyseal tibial joint block includes all the main fragments of the epiphysis and metaphysis.

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Variable angle-stable screw fixation in the metaphyseal part of the plate. Conventional or angle-stable screws placed in the diaphyseal part of the plate.

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
Longer plates on request to treat 2-level fractures.

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
When there is an intra-articular tibial fracture, the joint block is first restored and retained temporarily with wires that will not get in the way of the plate. Select a plate of suitable length for the appropriate side and position it strictly medially with the flange directed anteriorly until it fits optimally on the bone.

If necessary, shape the plate with bending instruments. The distal tab must always be bent also and if it is not needed it can be divided with a strong side cutter. The implant can be inserted minimally invasively. After screwing the drill sleeves into the most distal screw holes, the direction of the screws and the level of the plate's position relative to the joint can be set with a K-wire parallel to the drill sleeves. After a radiological check, fix the plate temporarily.

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The distal diverging screw holes are first filled with angle-stable screws by screwing in the drill sleeves, drilling with a 2.5 mm drill and inserting the small-fragment head-locking screws. Use the given screw directions as standard, particularly in the case of intra-articular fractures, in order to avoid collisions of the screws. The screws can be inserted in angle-stable mode with up to 20° deviation from the standard direction. This is usually necessary only with very small or large bones in order to ensure that the screw runs optimally from medial to lateral close to the cortex. Ensure that there are no collisions between the screws and consider using smaller cortical locking screws. The medial malleolus fragment is fixed with a conventional small-fragment screw. Check the outcome of the internal fixation critically by fluoroscopy. Irritation of the fibula and syndesmosis by screws inserted from the medial aspect must be avoided and pay special attention to the screw length.

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After reducing the restored joint block to the tibial shaft, the proximal plate holes are now filled with either angle-stable (as described above) or conventional small-fragment cortical screws. Use at least 7 cortical screws and 9 in the case of fractures extending far into the tibial shaft.

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
After the final radiological check and documentation, the wound is closed with a drain and a wellpadded lower leg splint until the wound has healed. Postoperative treatment is in accordance with the familiar guidelines. It is limited by the stability of the internal fixation and patient compliance.

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
Locating sieve (for instruments and implants)

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Hexagonal screw-driver, with handle, 2.5 mm, for screws diameter 3.5mm – 4.0 mm

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Drill guide for angle-stable systems, length 68 mm, SF

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Drill guide for variable angle-stable systems, length 68 mm, SF

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D 2,5 mm, Länge 155 mm

Twist drill (for quick coupler) diameter 2.5 mm, length 155 mm

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Kirschner wire with trocar point and round end, diameter 1.8 mm, length 150 mm

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*No part of this publication may be reproduced,
stored in a retrieval system or transmitted in any
form or by any means electronic,
mechanical, photocopying, recording or other-
wise without the prior written permission of the
publisher.*

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© Copyright Königsee Implantate GmbH
KÖNIGSEE IMPLANTATE GmbH