



medisun[®] Psori comb

EFFICIENT PHOTOTHERAPY OF THE SCALP
AND SMALL SKIN AREAS



Distributed in the UK by

ATHRODAX
HEALTHCARE INTERNATIONAL LTD

medisun® Psori comb

EFFICIENT PHOTOTHERAPY OF THE SCALP AND SMALL SKIN AREAS

medisun® Psori comb permits targeted treatment of small skin areas, e.g. on the hands, elbows, the hollows of the knees, etc. The comb row is used as a spacer and can be removed if necessary.

Extensive accessories included

The medisun® Psori comb is delivered in a padded transport box, with detailed usage instructions, UV protection goggles and a countdown timer.

The enclosed dosage-time table permits a reproducible dosage per second. The external shock-proof plug leads to a very light-weight and easy manageable lamp module for treating the scalp.



The padded transport box is practical and comfortable.



Product video
www.medisun.de/en


TECHNICAL DATA

DIMENSIONS in cm	L27 x D5.0 x H5.2
WEIGHT	approx. 0,9 kg
MAINS CONNECTION	230 V / 50 Hz
UV LAMPS	1 x 9 W
CERTIFICATION	CE 0197
WEEE-REG. NO.	DE 91065138


Technical changes reserved.

LAMP TYPE

UV range	UV lamps	Item-No.
UVA	1 x 9 W	3810
UVB ₃₁₁	1 x 9 W	3757



UVA



UVB-311



SCHULZE & BÖHM GmbH
Hermülheimer Str. 10
D-50321 Brühl
Germany

Distributed in the UK by
ATHRODAX
HEALTHCARE INTERNATIONAL LTD
Hawthorn Business Park, Drybrook, Glos. GL17 9HP
Telephone: 01594 544440 Fax: 01594 545800
Email: sales@athrodax.co.uk www.athrodax.co.uk